

# KL-383 OVER BED TABLE HEIGHT ADJUSTABLE

## Features:

Archimed overbed tables especially produced for hospital usage as C shape leg which allows to store bedside cabinet inside the shape C.

The construction is steel, top is ABS plastic easy to assemble and very practical . It has 4 castors, 2 of them has brake  
The C shape legs allows maximum approach to the bed for patient comfort.

Height is adjustable by gas spring. Height range is 84-110 cm.

ABS foot table dimension is 65cm x 40 cm

Metal profiles of ove bed are optional:

- a- Round
- b- Square



Code	Colour	Description	Height range (mm)	Table top (mm)	Width (mm)	Depth (mm)
KL-383_A	White	Overbed table with height adjustment	840 - 1100	650 X 400	650	400



Arşimed Medikal Kozm. San. Tic Ltd. Şti.

Mimar Sinan Mah. Bosna Cd. No:5 / 3 Çekmeköy - İstanbul

+90 216 567 8 567 (pbx) +90 216 567 8 367 (pbx)

www.arsimed.com.tr info@arsimed.com.tr

